



<b>Modular Arm Support Order Form</b>	Date:
Account number:	Phone:
P.O. number:	Email:
Company Name:	Clinician/Sales Rep/ATP:
Purchasing contact:	Marked for:

Dynamic Arm Support Kit		Part Number
Dynamic QR Arm Support, Small Hand Block Flex, Left		AS201L
Dynamic QR Arm Support, Small Hand Block Flex, Right		AS201R
Dynamic QR Arm Support, Large Hand Block Flex, Left		AS202L
Dynamic QR Arm Support, Large Hand Block Flex, Right		AS202R
Mounting Hardware (REQUIRED)		Part Number
Armrest Mounting Clamp, Ø19mm (Ø3/4")		AS119
Armrest Mounting Clamp, Ø22mm & 25mm (Ø7/8" & 1")		AS118
Armrest Mounting Plate		AS120
Modular Arm Support Components		Part Number
Dynamic QR Arm Support Base, Left		AS200L
Dynamic QR Arm Support Base, Right		AS200R
Small Hand Block Flex, Left		AS115L
Small Hand Block Flex, Right		AS115R
Large Hand Block Flex, Left		AS116L
Large Hand Block Flex, Right		AS116R
Assessment Kit		Part Number
Assessment Kit Includes Dynamic QR Arm Support Base, Small Hand Blocks, Large Hand Blocks, all Armrest Mounting Clamps and Armrest Mounting Plate.		AS001

### Hand Block Sizing Chart

	 <b>FINGER LENGTH</b> Measure length of middle finger from distal portion of the first phalangeal joint to tip of finger.	 <b>HAND WIDTH</b> Measure hand widths across widest part of first phalangeal joints.
Small	8.5cm - 10.5cm (3 3/8" - 4 1/8")	<= 9.0cm (3 1/2")
Large	10.5cm - 12.0cm (4 1/8" - 4 3/4")	<=11.0cm (4 1/4")